

About Our Company

HBW s.r.l. HUMAN BRAIN WAVE

- *HBW was founded in April 2012 to provide innovative medical services by producing biotechnological devices for medical applications.*
- *HBW was founded by Antonio Graziano and Riccardo D'Aquino following their research on stem cells and regenerative medicine.*
- *HBW in 2014 patented its micro-graft technology called RIGENERA[®].*
- *HBW started to test a prototype of Rigenera[®] technology in the Italian market in January 2013 (first clinical tests started in 2006)*
- *The Rigenera[®] technology is based on micrografts generated from autologous tissues without enzymes and/or additives using a disposable medical device called Rigeneracons/Adipecons.*
- *The Rigeneracons is CE and FMA (Japan) certified as a disposable medical device and FDA listed for Skin surgical graft and Bone mill*

ACSICON 2019
17th National Conference of
Association of Cutaneous Surgeons (I)

**11th WORLD CONGRESS
HAIR RESEARCH**
Sitges, Barcelona

AMWC

17th AESTHETIC & ANTI-AGING
MEDICINE WORLD CONGRESS

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
WORKSHOP
Rigenera
Re-evolution

to join us at
15, 2018
L - BANGKOK
1 Soi Ruamrudee
nwan, 10330

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
Rigenera worldwide: where we are where we'll go...


Rigenera® a HBW technology

The RIGENERACONS are surgical disposable medical devices CE, FMI approved and FDA listed in which the **autologous** tissue sample form the patient are **mechanically disaggregated** to obtain the Micrografts.

Every device owns a grid with 100 hexagonal holes. Any hole is embraced by six micro-blades designed for efficient cutting of hard and soft tissues.

How it works ?

- *After the tissue simple collection, RIGENERA allows a standard tissue disaggregation in one surgical time*

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Key points

- *In just one surgical time, the patient is the donor and the acceptor of calibrated micrografts of 80 Micron.*
- *The tissue is minimally manipulated since the devices only performs mechanical cutting and filtration.*
- *The collections of tissue particles in average size of 80 Microns, after filtration, allow the collection mainly of progenitor cells, extracellular matrix and growth factors.*
- *The micrografts can be deliverable in the acceptor site with a syringe and used alone or in association with different biomaterials according the tissue to regenerate.*

AUTOLOGOUS USE


HOMOLOGOUS TISSUE

ONE SURGICAL TIME

Rigenera® a HBW technology a step behind: What are the micrografts?


The micrografts idea is not new, Cicero Park Meek was pioneer of skin expansion techniques (1960) reaching a 10-fold skin expansion.

Proving that multiple pieces of a large graft provide more active edges for regeneration.

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- O conceito de utilização de microenxerto remonta à década de 1950(Cícero Parker Meek) Successful microdermagrafting using the Meek-Wall microdermatome. **Am J Surg. 1958 Oct;96(4):557-8.**

Rigenera® a HBW technology: The new Micrografts Theory

The Clinical Challenges:

To prepare viable micrografts in a fast, reliable, safe procedure without any chemical substance or tissue extensive manipulation

To reduce the size of the Micrografts particles and preserving as much as possible the cell viability

To overcome the Stem cells Legal, Ethical and Price limitations by using Micrografts-enriched in progenitor cells.

Editorial

Current Tissue Engineering, 2015, Vol. 4, No. 1 3

Editorial

New Perspectives and Therapies in Regenerative Medicine

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Tissue engineering is emerging as an interdisciplinary field in biomedical engineering that aims to repair diseased or damaged tissues or organs. The skin is the largest organ of the body in vertebrates and represents approximately one-tenth of the body mass in humans. It is composed of the epidermis and dermis with a complex nerve and blood supplying provides vital barrier function. Disruption of epidermal integrity due to trauma, disease, burn or surgery can be fatal, and therefore strategies are needed to enhance the physiological regenerative properties of the skin. Engineering skin substitutes represents a prospective source of advanced therapy in the clinical settings as

well as valuable skin surrogates for drug permeability tests and toxicity screening. However, at the present time, there are no models of bioengineered skin that completely replicate the anatomy, physiology, biological stability or aesthetic nature of uninjured skin.

Antonio Graziano

A promising alternative to tissue engineering is represented by stem cells therapy, in fact many studies have been conducted on its use to repair a damaged tissue or organ. The stem cell therapy is mainly based on the use of Mesenchymal Stem Cells (MSCs) which are multipotent adult stem cells with unique biological properties. Several in vitro studies and preclinical animal models reported that MSCs are promising for cell therapy showing the ability to home to sites of inflammation after tissue injury, to differentiate into various cell types and secrete multiple bioactive molecules capable of stimulating recovery of injured cells by inhibiting inflammation by a paracrine effect. MSCs also show the lack of immunogenicity and have the ability to exert immunomodulatory functions. Today, many pathological conditions are treated with MSCs, such as ischemic cardiovascular diseases, critical limb ischemia, bone and cartilage regeneration or neural diseases.

In this special issue, we aim to summarize some of the applications of MSCs, for example in facial aging, maxillo-facial defects and cell cardiomyoplasty. We also display a new approach to isolate specifically this cell population from human connective tissues and create autologous micro-grafts ready to clinical applications. Autologous bone grafting remains in fact a gold standard for the reconstruction for bone defects, mainly in the maxillofacial region, where dental pulp stem cells (DPSC) and bone-marrow mesenchymal stem cells (BMSC) represent the most common source of stem cells used for the building of 3D structures, due to their ability to self-renew and multi-lineage differentiation. In this context, was developed the Rigenera protocol, a new approach for regeneration of human injured tissues through the injection of autologous micro-grafts obtained by innovative medical device designed and produced for surgery. The patient is in fact the donor and acceptor of calibrated micro-grafts of 50 micron, displaying a high cell viability and an optimal regenerative potential. In the final part of this issue, we briefly show promising results obtained with Rigenera protocol mainly in the dentistry field.


Rigenera® a HBW technology: The science behind

- Principle 1: **The side population** .Numerous scientific papers have demonstrated that progenitor cells reside within a population with determined morphological features called side population. The main features of this side population are: 1) the size and 2) the expression of stem markers, significantly higher than in the wild population.
- Principle II : **The niche concept** : Preserving extra cellular matrix environment, the cells after desegregation still stand in their own physiological niche, that supports them in the grafting site.

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











What are made of the Rigenera Micrografts?

Rigenera Micro-grafts are 80/90 microns dimensions and they are made of


- High Viable cells
- Extra cellular matrix
- Growth Factors

Acceleration Mechanisms of Skin Wound Healing by
Autologous Micrograft in Mice
Shiro Jimi 1, *, Masahiko Kimura 2, Francesco De Francesco 3, Michele Riccio 3, Shuuji Hara 2
and Hiroyuki Ohjimi 4

Micrografts clinical application

AESTHETIC MEDICINE

- ◆ Androgenetic Alopecia
- ◆ Skin Rejuvenation
- ◆ Atrophic Scars
- ◆ Vitiligo
- ◆ Body and/or Face Remodeling with lipofilling

ODONTOLOGY

- ◆ Sinus Lift augmentation
- ◆ Socket Preservation
- ◆ Ulcers
- ◆ Maxillary bone necrosis
- ◆ Cleft lip and palate

WOUND HEALING

- ◆ Plastic and reconstructive surgery
- ◆ Diabetic foot
- ◆ Leg ulcers
- ◆ Non healing wound
- ◆ Burn wound
- ◆ Surgical dehiscence

ORTHOPEDICS /
REHABILITATION

- ◆ Cartilaginous defect
- ◆ Chondropathies
- ◆ Osteonecrosis
- ◆ Muscle Injury
- ◆ Tendon injury
- ◆ Ligament injury
- ◆ Meniscus
- ◆ Pseudarthrosis

Depending on the different application the tissue is collected from the patient by different anatomical region, in any cases the sample collection is an affordable and minimally invasive procedures

Scientific and Clinical Evidence Aesthetic Medicine DIVISION

HBW s.r.l. HUMAN BRAIN WAVE

An innovative regenerative treatment of scars with dermal micrografts

Fabiano Svolacchia, MD,^{1,*} Francesco De Francesco, MD,^{2,*} Letizia Trovato, Bsc, PhD,³ Antonio Graziano, BDsc, PhD,^{3,4,†} & Giuseppe A. Ferraro, MD, PhD^{2,†}

¹Centre of Nanoscience, Mavi Sud, Aprilia,

²Multidisciplinary Department of Medical-Surgical and Dental Specialties, Second University of Naples, Naples,

³Human Brain Wave srl, Turin, Italy

⁴SHRO – Temple University of Philadelphia, Philadelphia, PA, USA

Article

Platelet-Rich Plasma and Micrografts Enriched with Autologous Human Follicle Mesenchymal Stem Cells Improve Hair Re-Growth in Androgenetic Alopecia. Biomolecular Pathway Analysis and Clinical Evaluation

Pietro Gentile ^{1,*}, Maria G. Scioli ², Alessandra Bielli ², Barbara De Angelis ¹, Ciro De Sio ³, Domenico De Fazio ³, Gabriele Ceccarelli ^{4,5}, Angelo Trivisonno ³, Augusto Orlandi ², Valerio Cervelli ¹ and Simone Garcovich ⁶

Adipose Derived Stem Cells and Growth Factors Applied on Hair Transplantation. Follow-Up of Clinical Outcome

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Medical Science

Volume : 3 | Issue : 1 | Jan 2018 | e-ISSN : 2456-5040

CLINICAL AND HISTOLOGICAL EVALUATION OF THE RIGENERA[®] METHOD FOR THE TREATMENT OF ANDROGENETIC ALOPECIA

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ABSTRACT

Androgenetic alopecia has become a more common condition in society, affecting both genders. It is a disorder of multifactorial origin, with therapeutic options both in the rise and under development. Known options include the procedures of regenerative medicine with promising results. This paper assesses clinical and histological changes in patients with AGA after applying an autologous cellular suspension obtained using the Rigenera[®] system. After applying treatment, an increase in the mean of hair thickness, together with reduction of its loss, have been objectified; the level of satisfaction described by patients is worth noting. Based on the results, the improvement of AGA obtained with the Rigenera system is objective; these results need to be completed with data from future studies after using this promising technique.

Scientific and Clinical Evidence **Dental DIVISION**

HBW s.r.l. HUMAN BRAIN WAVE

Original Article

Sinus lift tissue engineering using autologous pulp micro-grafts: A case report of bone density evaluation

Giorgio Brunelli, Alessandro Motroni,¹ Antonio Graziano,¹ Riccardo D'Aquino,¹ Ilaria Zollino,¹ Francesco Carinci¹

Testo

Abstract:

Background: Although autografts are the standard procedure for bone grafting, the use of bone regeneration by means of dental pulp stem cell is an alternative that opens a new era in this field. Rigenera Protocol is a new technique able to provide the surgeon autologous pulp micro-grafts. **Materials and Methods:** At the Department of Oral Surgery, Don Orione Hospital, Bergamo, Italy, one patient underwent to sinus lift elevation with pulp stem micro-grafts gently poured onto collagen sponge. A CT scan control was performed after 4 month and DICOM data were processed with medical imaging software which gives the possibility to use a virtual probe to extract the bone density. Pearson's Chi-square test was used to investigate difference in bone density (BD) between native and newly formed bone. **Results:** BD in newly formed bone is about the double of native bone. **Conclusion:** This report demonstrated that micro-grafts derived from dental pulp poured onto collagen sponge are a useful method for bone regeneration in atrophic maxilla.

Key words:

Bone, homograft, jaw, reconstruction, resorption, stem cell

Dental Clinic, Don Orione Hospital, Bergamo, ¹Department of Morphology, Surgery and Experimental Medicine, University of Ferrara, Ferrara, Italy

Journal of Translational Science

Research Article

ISSN: 2059-268X

Periosteum-derived micro-grafts for tissue regeneration of human maxillary bone

Riccardo d'Aquino^{1,2}, Letizia Trovato^{2*}, Antonio Graziano², Gabriele Ceccarelli³, Gabriella Cusella de Angelis³, Angelo Marangini⁴, Alessandro Nisio⁵, Massimo Galli⁶, Massimo Pasi⁷, Marco Finotti⁸, Saturnino M Lupi⁹, Silvana Rizzo⁹ and Ruggero Rodriguez Y Baena⁹

¹Oral Surgeon - C.so Galileo Ferraris 63 Turin and Via R. G. De Ayala 9, Naples Italy

²Human Brain Wave srl, corso Galileo Ferraris 6, 10128 Turin, Italy

³Department of Public Health, Experimental Medicine and Forensics, University of Pavia, Pavia, Italy

⁴Oral Surgeon, Naples, Italy

⁵Oral Surgeon, Bari, Italy

⁶Oral Surgeon, Pistoia Italy

⁷Dentistry Service, San Raffaele Hospital, Milan, Italy

⁸Oral Surgeon, Padova, Italy

⁹School of Dentistry, Department of Clinico-Surgical, Diagnostic and Pediatric Sciences, University of Pavia, Pavia, Italy

Abstract

The bone regeneration is one of the most important challenges for regenerative medicine. In maxillo-facial area, bone resorption of the alveolar crest occurs after tooth extraction and leads to several risks for rehabilitation treatments, including dental implants procedures. Goal of our study was to demonstrate the efficacy of an innovative clinical protocol of bone tissue engineering called Rigenera protocol, aimed to create and optimize bio-complexes constituted by collagen biomaterial and human autologous periosteum-derived micro-grafts. We assessed the capacity of these bio-complexes to prevent the bone resorption in the alveolar crest with respect to simple collagen performing histological evaluations of neo-formed osseous tissue. We demonstrated that autologous bio-complexes significantly reduced the bone resorption of both horizontal and vertical dimension of alveolar crest when compared to collagen alone. We also showed that these bio-complexes accelerate the ossification process triggering the formation of new osseous tissue after 45 days from treatment and increasing the calcified matrix after 60 days and until to 120 days with respect to collagen alone. Taken together, these data showed the efficacy of bio-complexes composed by periosteum-derived micro-grafts and collagen in the alveolar ridge preservation through a reduction of bone resorption and an enhancement of new osseous tissue formation.

Scientific and Clinical Evidence Wound Healing DIVISION

HBW s.r.l. HUMAN BRAIN WAVE

Journal of Stem Cells Research, Reviews & Reports

Open Access

Austin Publishing Group

Case Report

A New Medical Device, Based on Rigenera Protocol, in the Management of Complex Wounds

Giaccone M^{1,2}, Brunetti M^{1,2}, Camandona M^{1,2}, Trovato L^{3*} and Graziano A³

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²Department of Surgical Sciences, University of Turin, Italy

³Human Brain Wave srl, Italy

*Corresponding author: Trovato L, Human Brain Wave, srl 10128 Turin, Italy, Email: info@hbwsrl.com, letizia.trovato@unito.it

Received: September 22, 2014; Accepted: October 22, 2014; Published: October 24, 2014

Abstract

Injury to the skin provides a unique challenge, as wound healing is a complex and intricate process. Acute wounds have the potential to move from into to chronic wounds, requiring the physician to have a thorough understanding of outside interventions to bring these wounds back into the healing cascade. The development of new and effective interventions in wound care remains an area of intense research. For this purpose, we aimed to report the use of a new medical device, called Rigeneracons® (Human Brain Wave srl, Turin), in the management of complex wounds of different etiology in two subjects where the wound healing process was very difficult. Rigeneracons® devices are based on Rigenera protocol that consists in the use of a cell population enriched of progenitor cells able to improve a tissue repair. By these case reports, we demonstrated the efficacy of Rigenera protocol in improving wound healing process through application of a cell suspension rich of progenitors cells in two different subjects that are donor and acceptor of these micro-grafts.

Keywords: Wound healing; Rigenera protocol; Rigeneracons®

jove Journal of Visualized Experiments

www.jove.com

Video Article

Tissue Characterization after a New Disaggregation Method for Skin Micro-Grafts Generation

Valeria Purpura¹, Elena Bondioli¹, Antonio Graziano², Letizia Trovato², Davide Melandri¹, Martina Ghetti¹, Andrea Marchesini³, Maria Gabriella Cusella De Angelis^{4,5}, Laura Benedetti^{4,5}, Gabriele Ceccarelli^{4,5}, Michele Riccio³

¹Burns Centre and Emilia Romagna Regional Skin Bank

²Human Brain Wave srl

³Plastic and Reconstructive Surgery, AO "Ospedali Riuniti"

⁴Department of Public Health, Experimental Medicine, Anatomy and Forensic, University of Pavia

⁵C.H.T Centre for Health Technologies, University of Pavia

Correspondence to: Antonio Graziano at lab@hbwsrl.com

URL: <http://www.jove.com/video/53579>

DOI: [doi:10.3791/53579](https://doi.org/10.3791/53579)

ANTICANCER RESEARCH 36: 975-980 (2016)

Treatment of Oncological Post-surgical Wound Dehiscence with Autologous Skin Micrografts

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¹Department of Plastic and Reconstructive Surgery, University of Turin,
San Giovanni Battista Hospital, Turin, Italy;

²Department of Medical Sciences, University of Turin, Turin, Italy;

³Unit of Orthopedics and Traumatology, Santa Croce Hospital, Moncalieri, Italy

Scientific and Clinical Evidence Orthopedic and Rehabilitation DIVISION

HBW s.r.l. HUMAN BRAIN WAVE

RESEARCH ARTICLE

Open Access

Rationale and pre-clinical evidences for the use of autologous cartilage micrografts in cartilage repair

Marco Viganò^{1†}, Irene Tessaro^{1†}, Letizia Trovato^{2*}, Alessandra Colombini¹, Marco Scala³, Alberto Magi⁴, Andrea Toto⁴, Giuseppe Peretti^{1,5} and Laura de Girolamo¹

OPEN

ORIGINAL ARTICLE

Reconstructive

Reconstruction of Alar Nasal Cartilage Defects Using a Tissue Engineering Technique Based on a Combined Use of Autologous Chondrocyte Micrografts and Platelet-rich Plasma: Preliminary Clinical and Instrumental Evaluation

pharmaceuticals

Article

In Vitro and In Vivo Studies of Alar-Nasal Cartilage Using Autologous Micro-Grafts: The Use of the Rigenera[®] Protocol in the Treatment of an Osteochondral Lesion of the Nose

Gabriele Ceccarelli^{1,2,*}, Pietro Gentile^{3,4}, Marco Marcarelli⁵, Martina Balli^{1,2}, Flavio Lorenzo Ronzoni^{1,2}, Laura Benedetti^{1,2} and Maria Gabriella Cusella De Angelis^{1,2,*}

Gentile et al., J Regen Med 2016, 5:2
DOI: 10.4172/2325-9620.1000129

Journal of
Regenerative Medicine

Research

A SCITECHNOL JOURNAL

A combined use of Chondrocytes Micro Grafts (CMG) Mixed with Platelet Rich Plasma (PRP) in Patients Affected by Pinch Nose Deformity


Gentile P^{1*}, Scioi MG², Bielli A³, Orlandi A³ and Cervelli V⁴

medium induced proliferation and a robust chondrogenesis of adult stem cells *in vitro* [5,6]. Direct mixing of chondrocytes with PRP leads to shrinking and deformed cartilage formation *in vivo* [7] owing to poor mechanical stability and rapid degradability. In particular the use of chondrocytes micrografts represents a micro-invasive procedure and grafts are more flexible to fill the lesions with various shapes. Today the main problem in transferring the experimental protocols of tissue engineering in the routine clinical practice is the identification of accessible sites where an adequate amount of stem cells are collected [8,9]. In addition, the need to specifically define technical procedure and its safety is an essential factor, as occurring for stem cell application in breast reconstruction and soft tissue defects [10-13].

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HISTOLOGICAL ANALYSIS SAFRANIN-O/FAST GREEN STAINING

CTR: Control
TR : Treatment with micrograft obtained from the knee joint

Advantages

1. Quick and easy application method that lasts 15 min
2. No hospitalization or operating room needed
3. Can be done in a single session
4. Only local anesthesia is necessary for punch extraction
5. Closed handling system
6. No extra equipment needed
7. No nursing or specialized assistant needed
8. Rapid recovery of the injury
9. Without side effects (Autologous Tissue)
10. Affordable price for the results obtained
11. No direct competitors

- Rigenera Micrografts for the management Androgenetic alopecia: Clinical - Italy. **ON GOING**
- Rigenera Micrografts for VITILIGO Repigmentation: Clinical Study - Italy. **PUBLISHED**
- Rigenera Micrografts for Chronic Ulcer : Double blind Case-Control Clinical Trial - Israel. **ON GOING – DATA ANALYSIS**
- Rigenera Micrografts for VITILIGO Repigmentation : Double blind Case-Control Clinical Trial - Israel. **ON GOING**
- Rigenera Micrografts for Surgical Dehiscence Prevention : Clinical Study - Italy. **ON GOING**
- Rigenera Adipose Micrografts for Enhancing the Lipofilling implantation - Italy . **COMPLETE – DATA ANALISYS**
- Rigenera Cartilage Micrografts the Management of Knee Chondropathy - Italy. **COMPLETE – DATA ANALISYS**
- Rigenera Cartilage Micrografts the Management of Knee Chondropathy - Japan. **ON GOING**
- Rigenera Cartilage Micrografts the Management of Knee Chondropathy - Italy. **COMPLETE – DATA ANALISYS**
- Rigenera Atrial micrograft transplantation for ischemic heart failure: Clinical Trial - Finland. **PARTIALLY PUBLISHED**
- Rigenera Micrografts enhance the healing in patient suffering of wide burns: Clinical Study – South Africa. **COMPLETE – PAPER SUBMITTED**
- Rigenera periosteum micrograft support bone regeneration in Socket preservation and Sinus lifts procedure: Clinical Trial - Brazil. **ON GOING**

HBW is partner of national and international research's projects

And the technology has been used for more than 60 Scientific publications

Autologous atrial micrograft (AAM) transplantation for ischemic heart failure coordinated by Dr. Esko Kankuri, University of Helsinki.

European Space Agency: **Tissue Healing in Space: Techniques for promoting and monitoring tissue repair and regeneration**

Piemonte Region- Bando: PIATTAFORMA TECNOLOGICA "SALUTE E BENESSERE" Progetto **Digital technology For Lung Cancer Treatment**

Rigenera Emergency Kit for military and civilian emergencies. The HBW is part of RAWINTS Project for Rapid Wound Healing, financed by NATO Science for Peace and Security Project.

EUROSTARS PROJECT : Redefine Regenerative medicine with a point of care Tissue Technology - **GRANTED**



<http://www.jhltonline.org>

Intraoperative processing and epicardial transplantation of autologous atrial tissue for cardiac repair

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
- Autologous atrial micrografts during coronary artery bypass grafting (CABG). Micrografts are isolated from the atrial appendage(AA) tissue and are combined with fibrin gel and a-decellularized small intestinal submucosa (SIS)Matrix sheet.
- The micrograft transplant was placed on the anterior wall of the left ventricle, on the scar area identified in the preoperative magnetic resonance imaging. The sheet was sutured to the myocardial surface with minimal disturbance.
- Magnetic resonance imaging at 6 months showed markedly improved movement and ventricular wall thickening in the treated infarcted scar area.

Rigenera micro grafts are effective in reducing the myocardial scar after infarction

STUDY PROTOCOL

Open Access

Epicardial delivery of autologous atrial appendage micrografts during coronary artery bypass surgery—safety and feasibility study

Annu Nurmi^{1*}, Tuomo Nieminen^{1,2}, Tommi Pätälä¹, Miia Lampinen⁴, Miia L. Lehtinen¹, Sari Kivistö³, Miia Holmström³, Erika Willman⁶, Kari Teittinen¹, Mika Laine¹, Juha Sinisalo¹, Markku Kupari¹, Esko Kankuri⁴, Tatu Juvonen¹, Antti Vento¹, Ralli Suojanen⁶, Ari Haajula¹ and The AADC consortium








Fig. 1 **a** Preparing the AADC-sheet. (A) The atrial appendage tissue is processed with cell therapy tissue homogenizer (Rigenera-system). (B) The micrografts are secured to extracellular matrix sheet (Cormatrix®) by using a fibrin sealant (Tisseel™). (C) The AADC sheet is placed to the myocardium in the location of infarction scar (animal model). **b** Administration of therapy during CABG surgery. Figure reproduced from our article by Lampinen et al. (Current Gene Therapy, 2015)

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